AMENDMENTS TO THE SPECIFICATION

At page 1, please amend the title as follows:

ASSAY FOR PREDICTING CELL ACTIVITY THE CAPACITY OF A CELL POPULATION TO INDUCE BONE FORMATION

At page 5, please amend the paragraph starting at line 29 as follows:

It is also possible to detect ALP by allowing it to convert a substrate for the enzyme ALP and detecting formed reaction product. Suitable substratea substrates in this respect are paranitro phenyl phosphate and alpha-naphtol AS-B1 alpha-naphthol AS-BI phosphate. Hydrolysis of the latter leads to the formation of a highly insoluble naphtol naphthol that may be coupled to a suitable diazonium salt that [[it]] is preferably present. A suitable example of such a diazonium salt is fast blue RR. The coupling provokes a coulour color reaction that can be detected by the naked eye or by UV. The reaction product obtained after conversion of paranitro phenyl phosphate by ALP may be detected using Sigma 104R SIGMA 104® phosphatase substrate (Sigma-Aldrich, St. Louis, MO) and UV. When the substrate is para-nitro phenyl phosphate, the cells are preferably first subjected to lysis and sonification.

At page 10, please amend the paragraph starting at line 14 as follows:

Human bone marrow stromal cells were plated at a density of 5000 cell per em(superscript: 2) cm² in 6 well plates (n=3) and cultured with and without dexamethasone for 7 days. After washing with PBS they were subjected to lysis and sonification. As substrate for ALP, para-nitro phenyl phosphate (PNP) was used (10 mM PNP in 1 ml diethanol amine and 1

USSN 10/730,770 Express Mail Receipt No. EV 515643376 US mM magnesium chloride (MgCl₂.6H₂O) at pH 9.8). The reaction product was detected using Sigma 104R SIGMA 104[®] phosphatase substrate (52.6 mg) dissolved in 10 ml ALP buffer (10 ml 10M diethanol amine and 90 ml demi-water to which, after overnight incubation, 10.33 mg magnesium chloride (MgCl₂.6H₂O) was added). The pH was adjusted to 9.8 using 1N HCl. 100 μl of this substrate was added to 100 μl of the cell lysate for 15 minutes at 37°C. A yellow colour color change was observed in a plate reader at 405 nm. The results obtained were calibrated and adjusted for the amount of DNA in a known manner.

At page 12, please amend the paragraph starting at line 29 as follows:

Human bone marrow stromal cells were seeded at a density of 5000 cells per em(superscript: 2) cm² in a T75 flask and cultured for 3-7 days in media with and without 1.25(OH)D3 (vitamin D3). After trypsinization the cells were blocked against nonspecific binding. Cells were then resuspended in FIX® solution A (fixation medium) (Caltag laboratories Laboratories, Burlingame, CA), incubated for 15 minutes and washed. A mixture 1:1 was made of PERM® solution B (permeabilization medium) (Caltag laboratories Laboratories, Burlingame, CA) and the osteocalcin mouse anti-human antibody (Zymed) dissolved in the blocking buffer (1:100 dilution). The cells were incubated for 15 minutes with the mixture. After incubation and washing, antibody reactivity was detected by suspending the cells with blocking buffer containing [[Goat]] goat anti-mouse FITC conjugated F(ab')2 F(ab')2 fragment (DAKO) (1:100 dilution) and incubating the cells for 30 minutes in the dark. After washing the cells were resuspended in FACS flow/staining and analyzed by a FACS-Calibur FACSCalibur apparatus (Beeton Diekenson BD Biosciences Immunocytometry systems, Franklin Lakes NJ). For each event 10.000 events were collected.